


HumaMeter A1c

Benchmark for Point of Care Testing

- › Accurate, reliable and fast HbA1c measurement
- › Easiest handling
- › Boronate Affinity Quenching Technology for most reliable results

CLINICAL CHEMISTRY



Human

Diagnostics Worldwide

HumaMeter A1c

Benchmark for Point of Care Testing

Boronate Affinity Method

- > Interference free
- > Results in only 4 minutes
- > Measuring range: 4–17% HbA1c
- > Imprecision: <3%
- > Sample volume: 4 µl!
- > Reagents stability for 30 days at RT

Easy handling


- > No pipetting capillary technology
- > Use of whole blood from finger-tip or venous blood samples
- > Innovative blood collector
- > Special designed cartridges
- > Step by step instructions on display

Out of routine features


- > Reporting in IFCC (mmol/mol) and NGSP/DCCT (%) values
- > 7000 results memory
- > Barcode reader for calibration code and patient data

HumaMeter A1c Technology

Boronate Affinity Quenching Technology is based on affinity of boronate for glycosylated proteins, fluorescent binding and fluorescent quenching.



- The HumaMeter A1c reagent has a fluorescent signal
- The reagent fluoresces strongly when unattached
- Fluorescence quenches when attached to HbA1c
- The fluorescence quenching observed is related to the quantity of HbA1c in the sample
- HumaMeter A1c can express the results in both IFCC and DCCT values



Ordering information

	REF
HumaMeter A1c	16080
Printer HumaMeter A1c	16081
HbA1c reagent kit 50 Tests	16085/50
HbA1c control	16086

References

- Concensus Committee, Consensus statement on the worldwide standardisation of the HbA1c measurement - Diabetologia, DOI 10.1007/s00125-007-0789-7
- Koenig BS, Peterson CM, Kilo C, Cerami A, Williamson JR - (1976) Hemoglobin A1C as an indicator of the degree of glucose intolerance in diabetes. Diabetes 25:230–232
- M. Silink et al, diabetes Voice 2007 ,12| Volume 52 | Issue 4